

RIH – CT ANGIOGRAM ABDOMEN/PELVIS S/P REPAIR GRAFT SIEMENS DEFINITION AS+ PROTOCOL

Indications: Evaluate patency of stent graft, to determine thrombosis of excluded portion of aorta, and to look for endovascular leaks.

Position/Landmark	Head first or feet first-Supine Sternal Notch																									
Topogram Direction	Craniocaudal / Craniocaudal																									
Respiratory Phase	Inspiration																									
Scan Type	Helical																									
Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization	Care kV 120 / Care Dose4D 180 / 0.5 sec 1.2:1 , 32.00mm 3 / 8																									
Detector width x Rows = Beam Collimation	0.625mm x 64 = 40mm (128 x .6mm)																									
Average Tube Output	ctdi – 10.0mGy dlp – 500mGy.cm																									
First Helical Set Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination .</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>nc abdomen/pelvis</td> <td>5mm x 5mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>2</td> <td>coronal nc abd/pelvis</td> <td>5mm x 5mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>3</td> <td>thin abd/pelvis</td> <td>.75mm x .6mm</td> <td>I40f medium</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination .	1	nc abdomen/pelvis	5mm x 5mm	I40f medium	pac	2	coronal nc abd/pelvis	5mm x 5mm	I40f medium	pac	3	thin abd/pelvis	.75mm x .6mm	I40f medium	terarecon					
recon	body part	thickness/ spacing	algorithm	recon destination .																						
1	nc abdomen/pelvis	5mm x 5mm	I40f medium	pac																						
2	coronal nc abd/pelvis	5mm x 5mm	I40f medium	pac																						
3	thin abd/pelvis	.75mm x .6mm	I40f medium	terarecon																						
Second Helical Set Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination .</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>axial ct angio</td> <td>3mm x 3mm</td> <td>I26f medium smooth</td> <td>pac</td> </tr> <tr> <td>2</td> <td>coronal ct angio</td> <td>3mm x 3mm</td> <td>I26f medium smooth</td> <td>pac</td> </tr> <tr> <td>3</td> <td>sagittal ct angio</td> <td>3mm x 3mm</td> <td>I26f medium smooth</td> <td>pac</td> </tr> <tr> <td>4</td> <td>thin ct angio</td> <td>.75mm x .6mm</td> <td>I26f medium smooth</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination .	1	axial ct angio	3mm x 3mm	I26f medium smooth	pac	2	coronal ct angio	3mm x 3mm	I26f medium smooth	pac	3	sagittal ct angio	3mm x 3mm	I26f medium smooth	pac	4	thin ct angio	.75mm x .6mm	I26f medium smooth	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination .																						
1	axial ct angio	3mm x 3mm	I26f medium smooth	pac																						
2	coronal ct angio	3mm x 3mm	I26f medium smooth	pac																						
3	sagittal ct angio	3mm x 3mm	I26f medium smooth	pac																						
4	thin ct angio	.75mm x .6mm	I26f medium smooth	terarecon																						
Third Helical Set Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination .</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>axial ct angio</td> <td>3mm x 3mm</td> <td>I26f medium smooth</td> <td>pac</td> </tr> <tr> <td>2</td> <td>coronal ct angio</td> <td>3mm x 3mm</td> <td>I26f medium smooth</td> <td>pac</td> </tr> <tr> <td>3</td> <td>sagittal ct angio</td> <td>3mm x 3mm</td> <td>I26f medium smooth</td> <td>pac</td> </tr> <tr> <td>4</td> <td>thin ct angio</td> <td>.75mm x .6mm</td> <td>I26f medium smooth</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination .	1	axial ct angio	3mm x 3mm	I26f medium smooth	pac	2	coronal ct angio	3mm x 3mm	I26f medium smooth	pac	3	sagittal ct angio	3mm x 3mm	I26f medium smooth	pac	4	thin ct angio	.75mm x .6mm	I26f medium smooth	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination .																						
1	axial ct angio	3mm x 3mm	I26f medium smooth	pac																						
2	coronal ct angio	3mm x 3mm	I26f medium smooth	pac																						
3	sagittal ct angio	3mm x 3mm	I26f medium smooth	pac																						
4	thin ct angio	.75mm x .6mm	I26f medium smooth	terarecon																						
Scan Start / End Locations	1 cm superior to diaphragm lesser trochanters 38cm																									
DFOV	decrease appropriately																									
IV Contrast Volume / Type / Rate	100mL Iohexol (Omnipaque 350) 4mL/sec																									
Scan Delay	Bolus tracking at level of celiac artery																									
2D/3D Technique Used	Workstream 4D mpr of 3mm x 3mm sagittal and coronal ct angiogram series, auto-transferred to PACS.																									
Comments:	Comments: A non-contrast study is done first. Then the cta is done using a smart prep at the level of the celiac artery. Note: There is a second helical scan done 60 seconds after the cta to look for subtle leak.																									
Images required in PACS	Topograms, 3mm x 3mm axial ct angio abdomen pelvis, 3mm x 3mm coronal ct angio abdomen pelvis, 3mm x 3mm sagittal ct angio abdomen pelvis, Patient Protocol																									