

## RIH – CHEST THORACIC ANEURYSM REPAIR SIEMENS DEFINITION AS+ PROTOCOL

**Indications:** Evaluate patency of stent graft, to determine thrombosis of excluded portion of aorta, and to look for endovascular leaks.

<b>Position/Landmark</b>	Head first or feet first-Supine 2cm superior to shoulders																				
<b>Topogram Direction</b>	Craniocaudal / Craniocaudal																				
<b>Respiratory Phase</b>	Inspiration																				
<b>Scan Type</b>	Helical																				
<b>Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization</b>	Care kV 120 / Care Dose4D 150 / 0.5 sec 1:1 , 24.00mm 3 / non con 4 contrast 7																				
<b>Detector width x Rows = Beam Collimation</b>	0.625mm x 64 = 40mm (128 x .6mm)																				
<b>Average Tube Output</b>	Each Helical: ctdi – 9 mGy dlp – 350 mGy.cm																				
<b>Helical Set</b> Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td><b>non con chest</b></td> <td>5mm x 5mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>2</td> <td><b>lungs</b></td> <td>5mm x 5mm</td> <td>I70f very sharp</td> <td>pac</td> </tr> <tr> <td>3</td> <td>thin chest</td> <td>.75mm x .7mm</td> <td>I40f medium</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination	1	<b>non con chest</b>	5mm x 5mm	I40f medium	pac	2	<b>lungs</b>	5mm x 5mm	I70f very sharp	pac	3	thin chest	.75mm x .7mm	I40f medium	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination																	
1	<b>non con chest</b>	5mm x 5mm	I40f medium	pac																	
2	<b>lungs</b>	5mm x 5mm	I70f very sharp	pac																	
3	thin chest	.75mm x .7mm	I40f medium	terarecon																	
<b>Helical Set</b> Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td><b>chest cta</b></td> <td>2mm x 2mm</td> <td>B30f medium smooth</td> <td>pac</td> </tr> <tr> <td>2</td> <td><b>coronal chest cta</b></td> <td>5mm x 5mm</td> <td>B30f medium smooth</td> <td>pac</td> </tr> <tr> <td>3</td> <td>thin chest</td> <td>.75mm x .7mm</td> <td>B30f medium smooth</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination	1	<b>chest cta</b>	2mm x 2mm	B30f medium smooth	pac	2	<b>coronal chest cta</b>	5mm x 5mm	B30f medium smooth	pac	3	thin chest	.75mm x .7mm	B30f medium smooth	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination																	
1	<b>chest cta</b>	2mm x 2mm	B30f medium smooth	pac																	
2	<b>coronal chest cta</b>	5mm x 5mm	B30f medium smooth	pac																	
3	thin chest	.75mm x .7mm	B30f medium smooth	terarecon																	
<b>Helical Set</b> Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td><b>delayed chest</b></td> <td>2mm x 2mm</td> <td>B30f medium smooth</td> <td>pac</td> </tr> <tr> <td>2</td> <td><b>coronal delayed</b></td> <td>5mm x 5mm</td> <td>B30f medium smooth</td> <td>pac</td> </tr> <tr> <td>3</td> <td>thin chest</td> <td>.75mm x .7mm</td> <td>B30f medium smooth</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination	1	<b>delayed chest</b>	2mm x 2mm	B30f medium smooth	pac	2	<b>coronal delayed</b>	5mm x 5mm	B30f medium smooth	pac	3	thin chest	.75mm x .7mm	B30f medium smooth	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination																	
1	<b>delayed chest</b>	2mm x 2mm	B30f medium smooth	pac																	
2	<b>coronal delayed</b>	5mm x 5mm	B30f medium smooth	pac																	
3	thin chest	.75mm x .7mm	B30f medium smooth	terarecon																	
<b>Scan Start / End Locations</b>  <b>DFOV</b>	1cm superior to lung apices mid kidney 38cm decrease appropriately																				
<b>IV Contrast Volume / Type / Rate</b>	100mL Iohexol (Omnipaque 350) / 4mL per second																				
<b>Scan Delay</b>	Bolus Tracking at the aortic arch																				
<b>2D/3D Technique Used</b>	Workstream 4D mpr of 5mm x 5mm <b>coronal chest</b> series of the arterial and delayed phases mip mode, auto-transferred to PACS.																				
<b>Comments:</b> Recon 3 is a thin helical volume of the chest that is archived to the TeraRecon server.																					
<b>Images required in PACS</b>	Topograms, 5mm x 5mm axial non con chest, 5mm x 5mm axial non con lungs, 2mm x 2mm axial arterial chest, 5mm x 5mm coronal arterial chest, 2mm x 2mm axial delayed chest, 5mm x 5mm coronal delayed chest, Patient Protocol																				