

**RIH – PEDI IV CONTRAST ABDOMEN/PELVIS AND LUMBAR SPINE
SIEMENS DEFINITION AS20 PROTOCOL**

Position/Landmark	Head first or feet first-Supine Sternal Notch																																			
Topogram Direction	Craniocaudal / Craniocaudal																																			
Respiratory Phase	Inspiration																																			
Scan Type	Helical																																			
Ref kV/Ref mAs/Rotation time (sec) Pitch / Speed (mm/rotation) Safire Strength / Dose Optimization	Care kV 120 / Care Dose4D 210 / 0.5sec .8:1 , 16.00mm 3 / 6																																			
Detector width x Rows = Beam Collimation	1.25mm x 16 = 20mm																																			
Average Tube Output	ctdi – 10.0mGy dlp – 500mGy.cm																																			
Helical Set Slice Thickness/ Spacing Algorithm Recon Destination	<table border="1"> <thead> <tr> <th>recon</th> <th>body part</th> <th>thickness/ spacing</th> <th>algorithm</th> <th>recon destination .</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>iv abdomen/pelvis</td> <td>3mm x 3mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>2</td> <td>coronal iv abd/pelvis</td> <td>3mm x 3mm</td> <td>I40f medium</td> <td>pac</td> </tr> <tr> <td>3</td> <td>axial lumbar spine</td> <td>3mm x 3mm</td> <td>I70f very sharp</td> <td>pac</td> </tr> <tr> <td>4</td> <td>coronal lumbar spine</td> <td>3mm x 3mm</td> <td>I70f very sharp</td> <td>pac</td> </tr> <tr> <td>5</td> <td>sagittal lumbar spine</td> <td>3mm x 3mm</td> <td>I70f very sharp</td> <td>pac</td> </tr> <tr> <td>6</td> <td>thin abd/pelvis</td> <td>1.5mm x 1mm</td> <td>I40f medium</td> <td>terarecon</td> </tr> </tbody> </table>	recon	body part	thickness/ spacing	algorithm	recon destination .	1	iv abdomen/pelvis	3mm x 3mm	I40f medium	pac	2	coronal iv abd/pelvis	3mm x 3mm	I40f medium	pac	3	axial lumbar spine	3mm x 3mm	I70f very sharp	pac	4	coronal lumbar spine	3mm x 3mm	I70f very sharp	pac	5	sagittal lumbar spine	3mm x 3mm	I70f very sharp	pac	6	thin abd/pelvis	1.5mm x 1mm	I40f medium	terarecon
recon	body part	thickness/ spacing	algorithm	recon destination .																																
1	iv abdomen/pelvis	3mm x 3mm	I40f medium	pac																																
2	coronal iv abd/pelvis	3mm x 3mm	I40f medium	pac																																
3	axial lumbar spine	3mm x 3mm	I70f very sharp	pac																																
4	coronal lumbar spine	3mm x 3mm	I70f very sharp	pac																																
5	sagittal lumbar spine	3mm x 3mm	I70f very sharp	pac																																
6	thin abd/pelvis	1.5mm x 1mm	I40f medium	terarecon																																
Scan Start / End Locations DFOV	1 cm superior to diaphragm lesser trochanters 38cm decrease appropriately																																			
IV Contrast Volume / Type / Rate	Contrast volume is 1cc per pound of body weight Omnipaque300 / 2cc per second or hand injection if necessary When oral contrast is prescribed, refer to the appropriate oral contrast agent's preparation and procedure guide.																																			
Scan Delay	50 seconds or just after hand bolus is completed																																			
2D/3D Technique Used	Workstream 4D mpr of 3mm x 3mm coronal abdomen/pelvis series, 3mm x 3mm axial lumbar series, 3mm x 3mm coronal lumbar series, 3mm x 3mm sagittal lumbar series, auto-transferred to PACS.																																			
Comments: Recon 6 is a thin helical volume of the abdomen/pelvis that is archived to the TeraRecon server.																																				
Images required in PACS	Topograms, 3mm x 3mm axial abdomen/pelvis, 3mm x 3mm coronal abdomen/pelvis, 3mm x 3mm axial lumbar spine, 3mm x 3mm coronal lumbar spine, 3mm x 3mm sagittal lumbar spine, Patient Protocol																																			